

Hydrogum

multiPro

Historical alginate for versatile use.



PRODUCT BENEFITS

- **High tear resistance:** helps to reduce the risk of the impression tearing when it is removed from the patient's mouth.
- **Elasticity:** optimum return to its original form



MINT

ONE OF OUR FIRST ALGINATES, ONE OF THE MOST POPULAR

Used in dental practices for over 20 years, Hydrogum was one of Zhermack's first alginates and is still one of the most popular in over 50 countries.

Clinical times**

Working time (including mixing time)	1' 10"
Time in mouth	1' 00"
Setting time	2' 10"

* Clinical use times are intended as from the start of mixing with deionised water at 23°C / 73°F.