

Zetalgin

easyPro

Essential alginate with a mint scent.



PRODUCT BENEFITS

- **Guaranteed quality:** meets ISO standards*, it has essential characteristics but affords guaranteed quality
- **Mint scent:** favours patient comfort



MINT

Clinical times**

Working time (including mixing time)	1' 35"
Time in mouth	1' 00"
Setting time	2' 35"

* ISO 21563:2013

** Clinical use times are intended as from the start of mixing with deionised water at 23°C / 73°F.

Zetalgin Chromatic

easyPro

Essential chromatic alginate with a mint scent.



PRODUCT BENEFITS

- **Guaranteed quality:** meets ISO standards*, it has essential characteristics but affords guaranteed quality.
- **Chromatic:** provides the practitioner with a visual guide during material mixing, working and positioning in the oral cavity.



MINT

The chromatic variation provides the practitioner with a visual guide during the different phases:

- purple during mixing
- grey during the working phase
- green during positioning in the patient's mouth

Clinical times**

Working time (including mixing time)	1' 10"
Time in mouth	1' 00"
Setting time	2' 10"

* ISO 21563:2013

** Clinical use times are intended as from the start of mixing with deionised water at 23°C / 73°F.